



Handheld 3-Function (AF-00300)

Specifications

Number of buttons	6
Number of functions	3, High/Low, leg spread and tilting function
Connector	Hooked connector 90°, 8-p DIN
Cable	Spiralcable 8x AWG26 1285mm long

Materials	PPC, TPE, SS
Housing colour	Gray
Dimension LxWxH	285x48x44mm

Order information

Article number	AF-00300
Replaces Article number	
Packaging	Per piece, 50pcs per box
M.O.Q.	10 pieces
Deliver time	3-5 working days



Suitable for patient hoist

+ Maxi Move

