



Charger 24V/DC SLA (AH-00804)

Specifications

Input voltage	100-240 VAC 50-60Hz
Input current	700-500mA
Output voltage	24 V/DC
Output current	400mA
IP class charger	IP40
Certification	UL2601-1 IEC60601-1 IEC60335-2-29
Operational temperature	0°C to +40°C
Storage temperature	-40°C to +70°C
Case colour	Black
Dimensions LxWxH	109x69x39mm

Order Information

Article number	C-AH-00801 for EU C-AH-00802 for GB C-AH-00803 for AUS C-AH-00804 for USA
Replacement for	KPA2000/KPA8200/SPL3024
Packaging	Individually packed, 50pcs per box
M.O.Q.	1 piece
Delivery Time	2-5 working days



Features

- + Suitable for AH-00200 /500 / KPA0100/ SPL3021 / KTA0102
- + Charger specifically designed for medical equipment
- + Delivered with power cord suited for UK, EU, USA or AUS
- + When the battery pack is charged, the charger will continue as trickle charger