



Charger 12V/DC SLA (C-AH-00904)

Specifications

Input voltage	100-240 VAC 50-60Hz
Input current	700-500mA
Output voltage	12 V/DC
Output current	500mA
IP class charger	IP40
Certification	UL2601-1 IEC60601-1 IEC60335-2-29
Operational temperature	0°C to +40°C
Storage temperature	-40°C to +70°C
Case colour	Black
Dimensions LxWxH	405x99x160mm



Order Information

Article number	C-AH-00901 for EU C-AH-00902 for GB C-AH-00903 for AUS C-AH-00904 for USA
Replacement for	NCA1000
Packaging	Individually packed, 50pcs per box
M.O.Q.	1 piece
Delivery Time	2-5 working days

Features

- + Suitable for Concerto / 102193
- + Charger specifically designed for medical equipment
- + Delivered with power adapter suited for UK, EU, USA or AUS
- + When the battery pack is charged, the charger will continue as trickle charger